



Signal processing

ERG



Group 17

Farida Mohammed : 24230035

Neama Refaat Hamed : 24230055

Table of Contents

01

Data Description

...

02

Roadmap

...

03


Data preprocessing


...

04

Data processing


...





01

Data Description



Patient's data

The data was gathered from both adults and pediatrics

Protocol	N. Adults	N. Pediatrics
Scotopic 2.0 ERG Response	23	53
Maximum 2.0 ERG Response	42	80
Photopic 2.0 ERG Response	32	74

The data consists of :

- The patient's age
- The Diagnosis
- The graph representing the response
- The potential difference after each half millisecond

...


A decorative graphic on the left side of the slide. It features a central blue circle with a white outline, connected by thin black lines to other blue circles of varying sizes. Some circles have white outlines, while others are solid blue. The lines and circles are arranged in a branching, network-like pattern. Ellipses (...) are placed near the top and bottom of the structure.

02


Roadmap


A decorative graphic on the right side of the slide, mirroring the one on the left. It consists of blue circles of different sizes connected by thin black lines in a branching pattern. Some circles have white outlines, and others are solid blue. Ellipses (...) are used to indicate the continuation of the structure.

Preprocessing roadmap




processing roadmap





03

Data Preprocessing



MR before organizing

Quantity	80 children, 42 adults					
Type	Maximum 2.0	Maximum 2.0	Maximum 2.0	Maximum 2.0	Maximum 2.0	Maximum
#	ERG Response	ERG Response	ERG Response	ERG Response	ERG Response	ERG Res
Diagnosis	368	360	360	363	363	
Age	The functional ac	The functional ac	The functional ac	The functional ac	The functional ac	The funct
	9.50	15.83	15.83	16.55	16.55	
	1 - 18	- children	18+	- adults		
a-wave amplitude, μV	41.40446617	36.36334441	39.29430409	41.73629335	43.85744978	44.45
b-wave amplitude, μV	52.30848728	60.34751823	74.18806126	55.66441905	53.13828863	41.1
a-wave latency, ms	21	21.5	21	20.5	21	
b-wave latency, ms	43.5	45	45	42.5	43	
a-wave amplitude, μV	28.63	60.29	- healthy children		41.61	50.1
b-wave amplitude, μV	55.69	83.19	- healthy children		89.68	104
a-wave latency, ms	13.29	28.71	- healthy children		16.86	17.4
b-wave latency, ms	39.08	56.92	- healthy children		37.94	39.9
Graph						
Time, ms	Signal, μV	Signal, μV	Signal, μV	Signal, μV	Signal, μV	Signal, μV
0	-3.335315569	-2.701023959	-4.064932455	7.181184174	1.120552042	25.42

A vertical line with several nodes. From top to bottom: a solid black dot, a solid black dot, a small blue circle with a white outline, three dots indicating continuation, a larger blue circle with a white outline, and a large blue circle with a white outline. A small solid black dot is located at the bottom left of the large blue circle.

The data frames were modified through:

- Removing null values
- Taking the important data only
- Mapping the diagnosis into normal and abnormal



...

Data after organizing

	Patient_id	Label	signal_type	img_Path
0	368	normal	MR	/kaggle/working/signals_MR/MR_1.png
1	360	normal	MR	/kaggle/working/signals_MR/MR_2.png
2	360	normal	MR	/kaggle/working/signals_MR/MR_3.png
3	363	normal	MR	/kaggle/working/signals_MR/MR_4.png
4	363	normal	MR	/kaggle/working/signals_MR/MR_5.png
...
267	362	normal	SR	/kaggle/working/signals_SR/SR_68.png
268	362	normal	SR	/kaggle/working/signals_SR/SR_69.png
269	365	abnormal	SR	/kaggle/working/signals_SR/SR_70.png
270	365	abnormal	SR	/kaggle/working/signals_SR/SR_71.png
271	366	abnormal	SR	/kaggle/working/signals_SR/SR_72.png


To get the signal:

- Plot the gathered data
- Convert the signals to images
- Apply Ricker wavelet transformation



Repetition and data augmentation:

- Repeated some of the signals
- Applied 5 methods for data augmentation
- Results: **1380 normal** and **1284 abnormal** images



Final dataframe

	img_Path	Label	img_id
0	/kaggle/working/augmented_images//MR_1_0_f_rot...	0	MR_1_0_f_rotated
1	/kaggle/working/augmented_images//MR_1_0_f_she...	0	MR_1_0_f_sheered
2	/kaggle/working/augmented_images//MR_1_0_f_zoo...	0	MR_1_0_f_zoomed
3	/kaggle/working/augmented_images//MR_1_0_f_hor...	0	MR_1_0_f_horizontal_flipped
4	/kaggle/working/augmented_images//MR_1_0_f_bri...	0	MR_1_0_f_brightness_changed
...
2683	/kaggle/working/signals/SR_66_1_f.png	0	SR_66_1_f
2684	/kaggle/working/signals/SR_68_1_f.png	0	SR_68_1_f
2685	/kaggle/working/signals/SR_69_1_f.png	0	SR_69_1_f
2686	/kaggle/working/signals/SR_72_1_f.png	1	SR_72_1_f
2687	/kaggle/working/signals/SR_73_1_f.png	1	SR_73_1_f



04

Data processing

Preparing data:

- Uploading the dataset
- Initializing the hyperparameters
⇒ image height & width, Batch size and n. Epochs
- Modifying the image paths
- Dividing data into train, test and validate datasets

train_df_aug

Label	Count
normal	1104
abnormal	1027

valid_df_aug

Label	Count
normal	138
abnormal	128

test_df_aug




Label	Count
normal	138
abnormal	129

Applying VGG16


- We modified the input layer so that we can specify the inputs
- For the outputs, we applied the sigmoid function to get binary results

...

Results of VGG16



Confusion matrix



Report

	precision	recall	f1-score	support
normal	0.55	0.62	0.59	138
retinal dystrophy	0.54	0.47	0.50	129
accuracy			0.55	267
macro avg	0.55	0.54	0.54	267
weighted avg	0.55	0.55	0.54	267

ROC curve

